


HINGED WRAPAROUND KNEE BRACE



Open Patella

Open patella and patella buttress help stabilize and protect kneecap



Triaxial Hinge

allows for unrestricted, natural knee movement



Wraparound Design

Sectional design and adjustable straps provide controlled compression and a custom fit

The Triaxial Hinge helps give you maximum medial-lateral support without restricting mobility.

The ultimate brace for weak or injured knees. The Triaxial Hinge closely simulates the natural tracking of the knee joint, while the open back design helps prevent bunching.

Recommended to help control subluxating patellas and to help relieve pain associated with Chondromalacia (irritated kneecap).

Support Level: maximum

Retail (clam)

3333REG Regular 3333LG Large 3333XL X-Large

Measure Around Knee

REG	LG	XL
12"-16"	16"-20"	20"-24"
30-40 cm	40-50 cm	50-60 cm